




SILICONE FOLEY CATHETER

Applications:

- Department of urology, anesthesia, ICU, department of general surgery, surgery, such as postoperative patients or patients that are unable to take care of themselves.



Features:

- It is made of imported medical grade silicone material.
- The silicone foley catheter has a larger inner lumen for better drainage than the same size one made from latex or PVC.
- No urate crystal and irritation occurs during intubation, thus catheter-associated urethral infection can be avoided.
- Silicone foley catheter is widely accepted because of better biocompatibility and indwelling period can be 30 days, which can reduce trauma to urethra caused by repeated intubation.

SILICONE FOLEY CATHETER

Order Information

Product Code	Description	Size(Fr)	Length(cm)	Balloon Capacity(cc)	Color
FX-020631	2 way	6	25	1.5-3mL	Green
FX-020831	2 way	8	31	3-5mL	Blue
FX-021031	2 way	10	31	3-5mL	Black
FX-021240	2 way	12	28	5-15mL	White
FX-021440	2 way	14	40	5-15mL	Green
FX-021640	2 way	16	40	5-30mL	Orange
FX-021840	2 way	18	40	5-30mL	Red
FX-022040	2 way	20	40	5-30mL	Yellow
FX-022240	2 way	22	40	5-30mL	Violet
FX-022440	2 way	24	40	5-30mL	Blue
FX-022640	2 way	26	40	5-30mL	Pink
FX-031640	3 way	16	40	5-30mL	Orange
FX-031840	3 way	18	40	5-30mL	Red
FX-032040	3 way	20	40	5-30mL	Yellow
FX-032240	3 way	22	40	5-30mL	Violet
FX-032440	3 way	24	40	5-30mL	Blue
FX-032640	3 way	26	40	5-30mL	Pink



GUANGDONG BAIHE MEDICAL TECHNOLOGY CO., LTD.

No. 89, Taoyuan East Road, Nanhai, Foshan, Guangdong Province, P.R.C.

Tel: +86-757-81207300 +86-757-81207302 +86-139-28979836

Fax: +86-757-81207311

Email: sales@baihemedical.com

Website: www.baihemedical.com



Version:202501