





## ® Central Venous Catheter Kit

---



***Focus on Quality, Respect for Life!***

# Central Venous Catheter



## Flexibility

Able Central Venous Catheter is made of superior polyurethane with excellent elasticity, flow rate is secured in venous obstructions.



## Kinking resistant

Guide wire is made of Nitinol with highly flexibility and kink resistance. It can be easily advanced even in venous obstruction. "J" and straight soft tip combined in one guide wire for your option.



















## Safety








Even and soft green tip reduces trauma to vessel, minimizing vessel erosion, hemothorax and cardiac tamponade.

# Completeness







Baihe provides full range of Central Venous Catheter in order to meet diversified clinical demands. Range from 24 G single lumen to 8.5 Fr Quad lumen.

Lumen	Needle	Hub	Guide wire
 Single			
 Double			
 Triple			
 Quad			



# Specification

Product Code	Specification	Guidewire	Dilator	Introducer Needle
FV-1121	Single lumen 24Ga 5cm	0.018"×50cm	4F×5cm	Straight/Y Needle 20Ga 
FV-11296	Single lumen 24Ga 6cm	0.018"×50cm	4F×5cm	
FV-1122	Single lumen 24Ga 8cm	0.018"×50cm	4F×5cm	
FV-1123	Single lumen 24Ga 10cm	0.018"×50cm	4F×5cm	
FV-1124	Single lumen 24Ga 13cm	0.018"×50cm	4F×5cm	
FV-1221	Single lumen 22Ga 5cm	0.018"×50cm	4F×5cm	Straight/Y Needle 20Ga 
FV-12296	Single lumen 22Ga 6cm	0.018"×50cm	4F×5cm	
FV-1222	Single lumen 22Ga 8cm	0.018"×50cm	4F×5cm	
FV-1223	Single lumen 22Ga 10cm	0.018"×50cm	4F×5cm	
FV-1224	Single lumen 22Ga 13cm	0.018"×50cm	4F×5cm	
FV-122915	Single lumen 22Ga 15cm	0.018"×50cm	4F×5cm	
FV-1225	Single lumen 22Ga 16cm	0.018"×50cm	4F×5cm	
FV-1226	Single lumen 22Ga 20cm	0.018"×50cm	4F×5cm	
FV-1321	Single lumen 20Ga 5cm	0.021"×50cm	4F×5cm	Straight/Y Needle 20Ga 
FV-13296	Single lumen 20Ga 6cm	0.021"×50cm	4F×5cm	
FV-1322	Single lumen 20Ga 8cm	0.021"×50cm	4F×5cm	
FV-1323	Single lumen 20Ga 10cm	0.021"×50cm	4F×5cm	
FV-1324	Single lumen 20Ga 13cm	0.021"×50cm	4F×5cm	
FV-132915	Single lumen 20Ga 15cm	0.021"×50cm	4F×5cm	
FV-1325	Single lumen 20Ga 16cm	0.021"×50cm	4F×5cm	
FV-1326	Single lumen 20Ga 20cm	0.021"×50cm	4F×5cm	
FV-1421	Single lumen 18Ga 5cm	0.025"×50cm	5.5F×6cm	Straight/Y Needle 19Ga 
FV-14296	Single lumen 18Ga 6cm	0.025"×50cm	5.5F×6cm	
FV-1422	Single lumen 18Ga 8cm	0.025"×50cm	5.5F×6cm	
FV-1423	Single lumen 18Ga 10cm	0.025"×50cm	5.5F×6cm	
FV-1424	Single lumen 18Ga 13cm	0.025"×50cm	5.5F×6cm	
FV-142915	Single lumen 18Ga 15cm	0.025"×50cm	5.5F×6cm	
FV-1425	Single lumen 18Ga 16cm	0.025"×50cm	5.5F×6cm	
FV-1426	Single lumen 18Ga 20cm	0.025"×50cm	5.5F×6cm	
FV-1428	Single lumen 18Ga 30cm	0.025"×70cm	5.5F×6cm	
FV-1523	Single lumen 16Ga 10cm	0.032"×50cm	6F×10cm	Straight/Y Needle 18Ga 
FV-1524	Single lumen 16Ga 13cm	0.032"×50cm	6F×10cm	
FV-152915	Single lumen 16Ga 15cm	0.032"×50cm	6F×10cm	
FV-1525	Single lumen 16Ga 16cm	0.032"×50cm	6F×10cm	
FV-1526	Single lumen 16Ga 20cm	0.032"×60cm	6F×10cm	
FV-1528	Single lumen 16Ga 30cm	0.032"×70cm	6F×10cm	
FV-1623	Single lumen 14Ga 10cm	0.032"×50cm	6F×10cm	Straight/Y Needle 18Ga 
FV-1624	Single lumen 14Ga 13cm	0.032"×50cm	8F×10cm	
FV-162915	Single lumen 14Ga 15cm	0.032"×50cm	6F×10cm	
FV-1625	Single lumen 14Ga 16cm	0.032"×50cm	8F×10cm	
FV-1626	Single lumen 14Ga 20cm	0.032"×60cm	8F×10cm	
FV-1628	Single lumen 14Ga 30cm	0.032"×70cm	8F×10cm	
FV/FC-2421	Double lumen 4F 5cm	0.018"×50cm	5.5F×6cm	Straight/Y Needle 20Ga 
FV/FC-24296	Double lumen 4F 6cm	0.018"×50cm	5.5F×6cm	
FV/FC-2422	Double lumen 4F 8cm	0.018"×50cm	5.5F×6cm	
FV/FC-2423	Double lumen 4F 10cm	0.018"×50cm	5.5F×6cm	
FV/FC-242912	Double lumen 4F 12cm	0.018"×50cm	5.5F×6cm	
FV/FC-2424	Double lumen 4F 13cm	0.018"×50cm	5.5F×6cm	

# Specification

Product Code	Specification	Guidewire	Dilator	Introducer Needle	
FV/FC-242915	Double lumen 4F 15cm	0.018"×50cm	5.5F×6cm	Straight/Y Needle 20Ga 	
FV/FC-2425	Double lumen 4F 16cm	0.018"×50cm	5.5F×6cm		
FV/FC-2426	Double lumen 4F 20cm	0.018"×50cm	5.5F×6cm		
FV/FC-2427	Double lumen 4F 25cm	0.018"×60cm	5.5F×6cm		
FV/FC-2428	Double lumen 4F 30cm	0.018"×70cm	5.5F×6cm		
FV/FC-2521	Double lumen 5F 5cm	0.021"×50cm	6F×7cm	Straight/Y Needle 20Ga 	
FV/FC-25296	Double lumen 5F 6cm	0.021"×50cm	6F×7cm		
FV/FC-2522	Double lumen 5F 8cm	0.021"×50cm	6F×7cm		
FV/FC-2523	Double lumen 5F 10cm	0.021"×50cm	6F×7cm		
FV/FC-2522912	Double lumen 5F 12cm	0.021"×50cm	6F×7cm		
FV/FC-2524	Double lumen 5F 13cm	0.021"×50cm	6F×7cm		
FV/FC-252915	Double lumen 5F 15cm	0.021"×50cm	6F×7cm		
FV/FC-2525	Double lumen 5F 16cm	0.021"×50cm	6F×7cm		
FV/FC-2526	Double lumen 5F 20cm	0.021"×50cm	6F×7cm		
FV/FC-2527	Double lumen 5F 25cm	0.021"×50cm	6F×7cm		
FV/FC-2528	Double lumen 5F 30cm	0.021"×50cm	6F×7cm		
FV/FC-2721	Double lumen 7F 5cm	0.035"×50cm	8F×10cm		Straight/Y Needle 18Ga 
FV/FC-2722	Double lumen 7F 8cm	0.035"×50cm	8F×10cm		
FV/FC-2723	Double lumen 7F 10cm	0.035"×50cm	8F×10cm		
FV/FC-2724	Double lumen 7F 13cm	0.035"×50cm	8F×10cm		
FV/FC-272915	Double lumen 7F 15cm	0.035"×50cm	8F×10cm		
FV/FC-2725	Double lumen 7F 16cm	0.035"×50cm	8F×10cm		
FV/FC-2726	Double lumen 7F 20cm	0.035"×60cm	8F×10cm		
FV/FC-2727	Double lumen 7F 25cm	0.035"×60cm	8F×10cm		
FV/FC-2728	Double lumen 7F 30cm	0.035"×70cm	8F×10cm		
FV/FC-272945	Double lumen 7F 45cm	0.035"×70cm	8F×10cm		
FV/FC-272960	Double lumen 7F 60cm	0.035"×80cm	8F×10cm		
FV-2823	Double lumen 8F 10cm	0.035"×50cm	9F×10cm	Straight/Y Needle 18Ga 	
FV/FC-2824	Double lumen 8F 13cm	0.035"×50cm	9F×10cm		
FV-292915	Double lumen 8.5F 15cm	0.035"×50cm	9F×10cm		
FV/FC-2825	Double lumen 8F 16cm	0.035"×50cm	9F×10cm		
FV/FC-2826	Double lumen 8F 20cm	0.035"×60cm	9F×10cm		
FV/FC-2828	Double lumen 8F 30cm	0.035"×70cm	9F×10cm		
FV-2924	Double lumen 8.5F 13cm	0.035"×50cm	9F×10cm	Straight/Y Needle 18Ga 	
FV-2925	Double lumen 8.5F 16cm	0.035"×50cm	9F×10cm		
FV-2926	Double lumen 8.5F 20cm	0.035"×60cm	9F×10cm		
FV-2928	Double lumen 8.5F 30cm	0.035"×70cm	9F×10cm		
FV/FC-3421	Triple lumen 4.5F 5cm	0.021"×50cm	6F×7cm	Straight/Y Needle 20Ga 	
FV/FC-342906	Triple lumen 4.5F 6cm	0.021"×50cm	6F×7cm		
FV/FC-3422	Triple lumen 4.5F 8cm	0.021"×50cm	6F×7cm		
FV/FC-3423	Triple lumen 4.5F 10cm	0.021"×50cm	6F×5cm 6F×7cm		
FV/FC-342912	Triple lumen 4.5F 12cm	0.021"×50cm	6F×5cm 6F×7cm		
FV/FC-3424	Triple lumen 4.5F 13cm	0.021"×50cm	6F×7cm		
FV/FC-342915	Triple lumen 4.5F 15cm	0.021"×50cm	6F×5cm 6F×7cm		
FV/FC-3425	Triple lumen 4.5F 16cm	0.021"×50cm	6F×7cm		
FV/FC-3426	Triple lumen 4.5F 20cm	0.021"×50cm	6F×7cm		
FV/FC-3427	Triple lumen 4.5F 25cm	0.021"×60cm	6F×5cm 6F×7cm		
FV/FC-3428	Triple lumen 4.5F 30cm	0.021"×70cm	6F×5cm 6F×7cm		

# Specification

Product Code	Specification	Guidewire	Dilator	Introducer Needle
FV/FC-3521	Triple lumen 5.5F 5cm	0.021"×50cm	6.5F×10cm	Straight/Y Needle 20Ga 
FV/FC-35296	Triple lumen 5.5F 6cm	0.021"×50cm	6.5F×10cm	
FV/FC-3522	Triple lumen 5.5F 8cm	0.021"×50cm	6.5F×10cm	
FV/FC-3523	Triple lumen 5.5F 10cm	0.021"×50cm	6.5F×10cm	
FV/FC-352912	Triple lumen 5.5F 12cm	0.021"×50cm	6.5F×10cm	
FV/FC-3524	Triple lumen 5.5F 13cm	0.021"×50cm	6.5F×10cm	
FV/FC-352915	Triple lumen 5.5F 15cm	0.021"×50cm	6.5F×10cm	
FV/FC-3525	Triple lumen 5.5F 16cm	0.021"×50cm	6.5F×10cm	
FV/FC-3526	Triple lumen 5.5F 20cm	0.021"×50cm	6.5F×10cm	
FV/FC-3527	Triple lumen 5.5F 25cm	0.021"×60cm	6.5F×10cm	
FV/FC-3528	Triple lumen 5.5F 30cm	0.021"×70cm	6.5F×10cm	
FV/FC-3721	Triple lumen 7F 5cm	0.035"×50cm	8F×10cm	
FV/FC-3722	Triple lumen 7F 8cm	0.035"×50cm	8F×10cm	
FV/FC-3723	Triple lumen 7F 10cm	0.035"×50cm	8F×10cm	
FV/FC-3724	Triple lumen 7F 13cm	0.035"×50cm	8F×10cm	
FV/FC-372915	Triple lumen 7F 15cm	0.035"×50cm	8F×10cm	
FV/FC-3725	Triple lumen 7F 16cm	0.035"×50cm	8F×10cm	
FV/FC-3726	Triple lumen 7F 20cm	0.035"×60cm	8F×10cm	
FV/FC-3727	Triple lumen 7F 25cm	0.035"×60cm	8F×10cm	
FV/FC-3728	Triple lumen 7F 30cm	0.035"×70cm	8F×10cm	
FV/FC-372945	Triple lumen 7F 45cm	0.035"×70cm	8F×10cm	
FV/FC-372960	Triple lumen 7F 60cm	0.035"×80cm	8F×10cm	
FV/FC-3923	Triple lumen 8.5F 10cm	0.035"×50cm	9F×10cm	Straight/Y Needle 20Ga 
FV/FC-3924	Triple lumen 8.5F 13cm	0.035"×50cm	9F×10cm	
FV/FC-392915	Triple lumen 8.5F 15cm	0.035"×50cm	9F×10cm	
FV/FC-3925	Triple lumen 8.5F 16cm	0.035"×50cm	9F×10cm	
FV/FC-3926	Triple lumen 8.5F 20cm	0.035"×60cm	9F×10cm	
FV/FC-3927	Triple lumen 8.5F 25cm	0.035"×60cm	9F×10cm	
FV/FC-3928	Triple lumen 8.5F 30cm	0.035"×70cm	9F×10cm	
FV-4923	Quad lumen 8.5F 10cm	0.035"×50cm	9F×10cm	
FV-492915	Quad lumen 8.5F 15cm	0.035"×50cm	9F×10cm	
FV-4924	Quad lumen 8.5F 13cm	0.035"×50cm	9F×10cm	
FV-4925	Quad lumen 8.5F 16cm	0.035"×50cm	9F×10cm	
FV-4926	Quad lumen 8.5F 20cm	0.035"×60cm	9F×10cm	
FV-4928	Quad lumen 8.5F 30cm	0.035"×70cm	9F×10cm	

## Remarks:

1. The product codes in this chart are for CVC kits with straight introducer needle configuration.
2. CVC kits with Y-shaped introducer needle configuration are also available in adult types. The product code of CVC kit with Y-shaped introducer needle configuration is to add letter "Y" after the product code of CVC kit with straight introducer needle configuration, for example, "FC-2726Y". Before you order, please clarify which configuration you request, either straight or Y-shape introducer needle.
3. Pediatric types of CVC include single lumen 22Ga, 20Ga, 18Ga; double lumen 4Fr, 5Fr; triple lumen 4.5Fr and triple lumen 5.5Fr.
4. Adult types of CVC include single lumen 16Ga, 14Ga; double lumen 7Fr, 8Fr, 8.5Fr; triple lumen 7Fr, 8.5Fr and quad lumen 8.5Fr.
5. "FV" in the product code for single lumen CVC means with cylinder hub.
6. "FC" in the product code for double lumen, triple lumen, quad lumen CVC means with triangular hub. "FV" in the product code for double lumen, triple lumen, quad lumen CVC means with round hub.



---

**GUANGDONG BAIHE MEDICAL TECHNOLOGY CO., LTD.**

No. 89, Taoyuan East Road, Nanhai, Foshan, Guangdong Province, P.R.C.

Tel: +86-757-81207300 +86-757-81207302 +86-139-28979836

Fax: +86-757-81207311

Email: [sales@baihemedical.com](mailto:sales@baihemedical.com)

Website: [www.baihemedical.com](http://www.baihemedical.com)



CMC Medical Devices & Drugs S.L.  
C/ Horacio Lengo No 18, CP 29006, Málaga, Spain.

